VarseoSmile TriniQ

Seite / page

de	Gebrauchsanweisung	2
en	Instructions for use	4
en-us	Instructions for use	6
fr	Notice d'utilisation	8
es	Instrucciones de uso	10
it	Istruzioni per l'uso	12
nl	Gebruiksaanwijzing	14
pt	Instruções de utilização	16
ru	Инструкция по применению	18
bg	Инструкция за употреба	20
cs	Návod k použití	22
da	Brugsanvisning	24

Seite / page

el	Οδηγίες χρήσης	26
et	Kasutusjuhend	28
hr	Upute za uporabu	30
hu	Használati utasítás	32
lt	Naudojimo instrukcija	34
lv	Lietošanas pamācība	36
pl	Instrukcja użytkowania	38
ro	Instrucțiuni de utilizare	40
sk	Návod na použitie	42
sl	Navodila za uporabo	44
sv	Bruksanvisning	46
tr	Kullanma talimatı	48

BEGO Bremer Goldschlägerei Wilh. Herbst GmbH & Co. KG
Wilhelm-Herbst-Str. 1 · 28359 Bremen, Germany
Tel. +49 421 20 28-0 · Fax +49 421 20 28-100
E-Mail: info@bego.com · www.bego.com

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VarseoSmile TriniQ

1. Produktbeschreibung

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VarseoSmile TriniQ

1. Device description

Ceramic-filled, acrylate-based hybrid material for the fabrication of permanent or temporary dental restorations for the treatment of tooth defects or tooth loss, including single tooth restorations (e.g. crowns, inlays, onlays, veneers), bridges, veneering on metal frames and artificial teeth. The hybrid material is moulded to create the restoration using the 3D printing process (SLA, DLP, LCD) and cured in a UV light-curing unit. The 3D printing and curing are carried out at light wavelengths of 385 nm to 405 nm.

2. Intended purpose

Resin for 3D printing of single tooth restorations, bridges and artificial teeth.

3. Indications and clinical benefit

The use of dental restorations made from VarseoSmile TriniQ is indicated in the following cases for restoring masticatory function and aesthetics of teeth in the anterior and posterior regions, including the occlusal surfaces:

- Tooth defects caused by trauma, caries, wear, abrasion or erosion/ biocorrosion
- Weakening and/or incomplete development of the tooth structure (genetic dysplasia)
- Occlusal changes, shape and positional changes
- Partial or complete edentulism

4. Users and target patient group

VarseoSmile TriniQ is a medical device that is intended for use by dentists, dental surgery assistants and dental technicians and can be used to treat all patients with no contraindications as defined in Section 5.

5. Contraindications

VarseoSmile TriniQ is contraindicated in the case of known allergies to one or more ingredients. If in doubt, a specific test should be done before using this device to confirm or rule out the particular allergy. The material should not be used for purposes other than those described by the manufacturer. Any deviation from the intended purpose, the instructions for use or the design requirements for the specific purpose can have negative effects on the chemical and physical quality of the objects fabricated using the material.

6. Material properties and pack size

VarseoSmile TriniQ is a ceramic-filled hybrid material made from (meth)acrylates, ceramic fillers, initiators, additives and pigments. VarseoSmile TriniQ is a polymer-based material that is cured by the application of UV light and contains photopolymerisation initiators for this purpose (Class 2 material in accordance with ISO 4049:2019).

Physical properties

Flexural strength	120 MPa
Surface hardness	>90 Shore D
Wavelength	385 – 405 nm

Shades and pack sizes

41170	VarseoSmile TriniQ	A1	500 g
41171	VarseoSmile TriniQ	A2	500 g
41172	VarseoSmile TriniQ	A3	500 g
41177	VarseoSmile TriniQ	A3.5	500 g
41173	VarseoSmile TriniQ	B1	500 g
41178	VarseoSmile TriniQ	B2	500 g
41174	VarseoSmile TriniQ	C2	500 g
41179	VarseoSmile TriniQ	D3	500 g
41175	VarseoSmile TriniQ	OM1	500 g
41176	VarseoSmile TriniQ	OM3	500 g
41180	VarseoSmile TriniQ	A1	250 g
41181	VarseoSmile TriniQ	A2	250 g
41182	VarseoSmile TriniQ	A3	250 g
41187	VarseoSmile TriniQ	A3.5	250 g
41183	VarseoSmile TriniQ	B1	250 g
41188	VarseoSmile TriniQ	B2	250 g
41184	VarseoSmile TriniQ	C2	250 g
41189	VarseoSmile TriniQ	D3	250 g
41185	VarseoSmile TriniQ	OM1	250 g
41186	VarseoSmile TriniQ	OM3	250 g

Shades and shade names based on VITA Classical A1-D4® shade guide and VITA Toothguide 3D-MASTER®.

Note: The availability of individual product versions may vary by region. Up-to-date information for each product can be found on the BEGO website.

7. Warnings

The following warnings must be observed when using VarseoSmile TriniQ in liquid form:

- Read the safety data sheet before use.
- VarseoSmile TriniQ must only be used for the fabrication of permanent or temporary bridges, crowns, inlays, onlays, veneers and artificial teeth following the design requirements described in these instructions for use.
- Any deviation from the instructions for use can have a negative effect on the chemical and physical qualities and the biocompatibility of the end product.
- The device contains substances that can cause skin irritation or allergic skin reactions.
- Avoid breathing in or swallowing. May lead to irritation of the respiratory passages.
- May harm fertility or the unborn child.
- If the resin comes into contact with the skin, wash with plenty of water.
- In the event of eye contact, rinse cautiously with water for several minutes. Remove contact lenses, if worn and if easy to do. Continue rinsing.
- If the product is swallowed, or if this is suspected, seek medical advice and call an emergency poison centre/a doctor if you feel unwell.

8. Precautions

Always wear protective clothing, safety goggles and nitrile gloves when handling VarseoSmile TriniQ. Further information on handling the product can be found in the safety data sheet, which is available at www.bego.com. We cannot completely rule out the possibility that reactions (e.g. intolerances or allergic reactions) could occur in isolated cases. If this occurs, the user should refrain from using VarseoSmile TriniQ.

Store in a well-ventilated location. Keep container closed tight. Keep locked away. Observe the storage and transport conditions. Do not use the resin after the expiration date.

These instructions for use or the relevant information contained they contain must be passed on to other users, as failure to comply with the instructions poses a high risk of fabricating an unreliable and/or unusable product that could put the user's safety at risk.

9. MR safety

VarseoSmile TriniQ is made from materials that are not electrically conductive, metallic or magnetic.

10. Adverse side-effects

Adverse side-effects, reactions or similar incidents, including any not listed in these instructions for use, must be immediately reported to the manufacturer or the local distributor.

11. Transport and storage

VarseoSmile TriniQ resin must be protected against light to prevent spontaneous polymerisation. It must be stored and transported in the original bottle at a temperature between 4°C and 35°C. Residual resin can be left in the resin tank of the 3D printer if the printer's UV light protection cover is kept closed or if the resin tank has a lightproof cover. The expiration date on the label of each original bottle specifies when the resin must be used by.

Note: If the material is used after the expiration date, the results may be unusable.

Fully cured printed objects can be stored at room temperature and protected against light.

12. Disposal

Fully cured material residue (base plate, support structure) can be disposed of with household waste. Unused resin or cleaning solutions containing resin residue must be disposed of according to local and national regulations.

13. How to use VarseoSmile TriniQ

VarseoSmile TriniQ has been developed for the fabrication of dental restorations using the 3D printing method. An overview of all 3D printers, cleaning agents and UV light-curing units that have been validated and approved for use with VarseoSmile TriniQ can be found at: <https://www.bego.com/de/3d-druck/kompatibilitaetsuebersicht/>

13.1. Required equipment

Warning: It is only possible to achieve results that reliably meet the intended purpose of VarseoSmile TriniQ by using 3D printers, cleaning agents and UV light-curing units that have been validated and approved for processing VarseoSmile TriniQ, and with the processing parameters specified for VarseoSmile TriniQ.

- Validated 3D printer, including nesting software, e.g. BEGO Varseo XS and BEGO CAM Creator
- Validated cleaning agent/device for the removal of residual resin, e.g. ethanol (≥96 %), isopropanol (≥99.5 %) or InovaPrint Wash* (REF 41167) and ultrasonic bath
- Validated UV light-curing unit, e.g. Varseo Cure
- VarseoSmile TriniQ
- Stainless steel spatula
- Spray bottle of ethanol or isopropanol
- Cutting disc or side cutters (to remove the support structure)
- Possibly sandblaster 1.5 bar
- Possibly blasting material 50 µm (e.g. Perlablast® micro, BEGO, REF 46092/54302)
- Protective equipment

13.2. Design requirements

The following information must be observed for the digital design of dental restorations made from VarseoSmile TriniQ:

Indication	Design	Image**	Permanent (12 months or more)	Temporary (up to 12 months)	Design requirement
Single tooth restorations	Fully anatomical single crown in the anterior region		✓	✓	<ul style="list-style-type: none"> • Minimum wall thickness for posterior and anterior teeth 0.7 mm • Minimum wall thickness for veneers 0.5 mm
	Fully anatomical single crowns in the posterior region		✓	✓	
	Veneer		✓	✓	
	Inlay, onlay		✓	✓	
Bridge constructions	Fully anatomical monolithic bridge constructions for the anterior and/or posterior region, up to 7 units without pontics		✓	✓	<ul style="list-style-type: none"> • Maximum width of the pontic 1 molar • Bridges in the anterior region: <ul style="list-style-type: none"> • Minimum wall thickness occlusal and circumferential 1.0 mm • Minimum connector cross section 14 mm² • Bridges in the posterior region: <ul style="list-style-type: none"> • Minimum wall thickness occlusal and circumferential 1.0 mm • Minimum connector cross section 16 mm²
	Fully anatomical monolithic bridge constructions for the anterior and/or posterior region with max. one pontic, up to 3 units		✓	✓	
	Fully anatomical monolithic bridge for the anterior and/or posterior region with max. 3 non-neighbouring pontics (up to 7-unit bridge constructions)		✗	✓	
Artificial teeth	Artificial teeth for connection with a denture base		✓	✓	<ul style="list-style-type: none"> • Minimum wall thickness 1.0 mm

** images shown are examples

13.3. Nesting

Required software, depending on the 3D printer being used:

- CAM Creator Print version 1.32 or higher
- Asiga Composer version 1.3.3 or higher

The printing parameters required for the processing of VarseoSmile TriniQ are given in the specific material file for your 3D printer. To select the correct material file, select your 3D printer model and VarseoSmile TriniQ in the corresponding nesting software. Load the desired STL file into your nesting software.

Additional printing parameters:

- Layer thickness 50 µm (preset/default)
- Horizontal alignment of the dental restorations, occlusal plane to the build platform. Exception for bridges: Align bridges at an angle of 45° to the build platform
- Automatic generation of support structures is recommended

* This symbol is a commercial designation/registered trademark of a company that is not part of the BEGO company group.

13.4. 3D Printing

The ideal processing temperature for VarseoSmile TriniQ is between 18°C and 35°C. Shake the bottle for approx. 2 minutes before first use to homogenise the resin. When filling the 3D printer's resin tank, make sure that the resin is only exposed to daylight as briefly as possible. Mix the resin in the resin tank if a transparent layer is visible on the surface.

Start the 3D print job following the instruction manual provided by the manufacturer of the 3D printer. Remove the printed objects from the build platform according to the instruction manual provided by the manufacturer of the 3D printer. Use a stainless steel spatula so if necessary.

Note: Once the printed objects have been removed from the 3D printer and separated from the build platform, the following steps, namely removal of residual resin and post-curing, should be carried out without any lengthy delay or storage of the printed objects so as to prevent undesirable material changes.

13.5. Removal of residual resin

Residual resin is removed by cleaning the printed objects in an ultrasonic bath. Ethanol (≥96%), isopropanol (≥99.5%) or InovaPrint Wash* (REF 41167) can be used as the cleaning solution:

1. 3 minutes cleaning in an ultrasonic bath in reusable cleaning solution.
2. 2 minutes cleaning in an ultrasonic bath in fresh cleaning solution.
3. Spray the printed objects with ethanol (≥96%) or isopropanol (≥99.5%).
4. Carefully dry the cleaned objects with compressed air.
5. If shiny spots remain on the printed objects, this indicates that there is still residual resin. Remove this by spraying the printed objects with ethanol (≥96%) or isopropanol (≥99.5%). Stubborn residues can also be removed with a brush soaked in ethanol (≥96%) or isopropanol (≥99.5%).
6. If necessary, dry the cleaned objects again with compressed air.

Warning: Never put ethanol or isopropanol into the ultrasonic bath itself. Always put it in the recommended container (REF 19621) and place the container in the ultrasonic bath filled with water. Use an explosion-proof ultrasonic bath.

Note: The specified total cleaning time should not be exceeded, as this could have a detrimental effect on the printed objects.

Note: When using ethanol (≥96%) or isopropanol (≥99.5%) as the cleaning solution, a whitish layer may form on the surface of the dental restorations. This is a deposit of the ceramic filling in VarseoSmile TriniQ. The deposit is harmless and can be removed by polishing the printed objects after post-curing.

13.6. Post-curing

Remove the support structures using a cutting disc or side cutters prior to post-curing. Take care not to deform the printed objects.

Post-cure the printed objects using the Varseo Cure UV light-curing unit with the parameters stored for VarseoSmile TriniQ for 2 x 2:30 min. Rotate the objects between the post-curing cycles to ensure even polymerisation.

An overview of other UV light-curing units that have been validated and approved for post-curing dental restorations made from VarseoSmile TriniQ can be found at:

<https://www.bego.com/de/3d-druck/kompatibilitaetsuebersicht/>

Note: The specified post-curing time only applies to UV light-curing units that are maintained on a regular basis.

Note: There may be a slight yellow discolouration of the printed objects after post-curing. This discolouration is normal and will soon disappear again.

Warning: If the post-curing process is interrupted due to a malfunction of the UV light-curing unit, the printed object must not be used until it has been fully cured with a complete post-curing cycle.

13.7. Finishing and further processing of printed objects

Post-cured objects made from VarseoSmile TriniQ can be polished with pumice stone and polishing compound. When polishing, take care not to deform the printed objects.

Post-cured objects made from VarseoSmile TriniQ can be customised or finished and supplemented with composite stains as well as direct and indirect composites. The current instruction manual from the manufacturer of the composite stains or composite veneering material must be observed.

Attaching crowns, inlays, onlays and veneers

- Permanent attachment of restorations made from VarseoSmile TriniQ onto a natural tooth can be done using fully adhesive or self-adhesive luting composites.
- Temporary attachment of restorations made from VarseoSmile TriniQ onto a natural tooth can be done using temporary luting cements.
- Attachment of restorations made from VarseoSmile TriniQ onto metal frames can be done using luting composites or veneering composites.

Note: The instructions of the relevant composite or cement manufacturer must always be followed.

Note: Glass ionomer cements are not suitable for attaching restorations made from VarseoSmile TriniQ.

Attachment of artificial teeth

Artificial teeth made from VarseoSmile TriniQ can be connected with cold-cured or 3D printed denture base plastics. Follow the current instructions from the manufacturer of the denture base material for how to connect the artificial teeth to the denture base. If the manufacturer of the denture base has not provided any instructions, the denture base material can be used to connect the artificial teeth to the denture base in the case of cold-cured denture bases. For 3D printed denture bases, the liquid VarseoSmile TriniQ resin can be used.

Cleaning post-cured objects

Fully post-cured objects made from VarseoSmile TriniQ can be steam-cleaned or immersed in a bath of cleaning solution.

14. Label symbols

Ref. number and symbol	Title	Description	Reference
	CE conformity mark	Indicates conformity with EU Regulation 2017/745	EU Regulation 2017/745 Appendix V
5.1.1 	Manufacturer	Indicates the manufacturer of the medical device	ISO 15223-1:2021 Medical devices • Symbols to be used in the information provided by the manufacturer
5.1.3 	Date of manufacture	Indicates the date when the medical device was manufactured	ISO 15223-1:2021 Medical devices • Symbols to be used in the information provided by the manufacturer
5.1.4 	Use-by date	Indicates the date after which the medical device must no longer be used	ISO 15223-1:2021 Medical devices • Symbols to be used in the information provided by the manufacturer
5.1.5 	Batch code	Indicates the manufacturer's batch code for identifying the batch or lot	ISO 15223-1:2021 Medical devices • Symbols to be used in the information provided by the manufacturer
5.1.6 	Catalogue number	Indicates the manufacturer's catalogue number for identifying the medical device	ISO 15223-1:2021 Medical devices • Symbols to be used in the information provided by the manufacturer
5.3.2 	Keep away from sunlight	Indicates a medical device that needs to be protected against sources of light	ISO 15223-1:2021 Medical devices • Symbols to be used in the information provided by the manufacturer
5.3.7 	Temperature limits	Indicates the temperature limits to which the medical device can safely be exposed	ISO 15223-1:2021 Medical devices • Symbols to be used in the information provided by the manufacturer
5.4.3 	Follow the printed or electronic instructions for use	Informs the user that it is necessary to observe the instructions for use	ISO 15223-1:2021 Medical devices • Symbols to be used in the information provided by the manufacturer
5.4.4 	Caution	Indicates that care must be taken when handling the device or the control unit in the area near where the symbol is shown, or that the current situation requires caution on the part of the operator or an intervention by the operator in order to prevent unwanted consequences	ISO 15223-1:2021 Medical devices • Symbols to be used in the information provided by the manufacturer
5.7.7 	Medical device	Indicates that the object in question is a medical device	ISO 15223-1:2021 Medical devices • Symbols to be used in the information provided by the manufacturer

* This symbol is a commercial designation/registered trademark of a company that is not part of the BEGO company group.

Instructions for use

VarseoSmile TriniQ

Caution: Federal law restricts this device to sale by or on the order of a dentist. For professional use only – Rx only.

1. Product description

Ceramic-filled, acrylate-based hybrid material for the fabrication of dental objects for the definitive or temporary restoration of dental defects or tooth loss in the form of single-tooth restorations (e.g. crowns, inlays, onlays, veneers), bridges and denture teeth. The hybrid material is formed into the respective dental object using the 3D printing process (SLA, DLP, LCD) and cured in a UV polymerization unit. The 3D printing as well as the polymerization take place at light wavelengths of 385 nm to 405 nm.

2. Purpose

Resin for 3D printing of single tooth restorations, bridges and artificial teeth.

3. Indications for use

VarseoSmile TriniQ is indicated as an indirect restorative for both anterior and posterior restorations, including occlusal surfaces. The VarseoSmile TriniQ material is used for fabricating permanent restorations such as inlays, onlays, veneers, full crown and bridge restorations. VarseoSmile TriniQ can also be used for the fabrication of artificial teeth and temporary crowns & bridges.

4. Users and target patient group

VarseoSmile TriniQ is a medical device intended for use by professionals, such as dentists, dental assistants and dental technicians, and can be used for the treatment of all patients who do not have a contraindication according to paragraph 5.

5. Contraindications

VarseoSmile TriniQ is contraindicated in case of known allergies to one or more ingredients. In cases of doubt, the allergy should be clarified and excluded by a specific test before using this product. The material should not be used for purposes other than those described by the manufacturer. Any deviation from the intended use, instructions for use or design requirements for the specific purpose may have a negative impact on the chemical and physical quality of the items made from this material.

6. Material properties and supply units

VarseoSmile TriniQ is a ceramic-filled hybrid material consisting of (meth-)acrylates, ceramic fillers, initiators, additives, and pigments. VarseoSmile TriniQ contains VarseoSmile TriniQ is a polymer-based material whose setting is achieved by the extra-oral application of UV-light and contains photopolymerization initiators for this purpose (Type 1, Class II, Group II material according to ISO 4049:2019).

Material properties

Flexural strength	120 MPa
Surface hardness	>90 Shore D
Wavelength	385 – 405 nm
Curing time (post curing)	See section Post Curing
Intensity for curing	>5 mW/cm ²
Depth of cure	>99 %

Supply units

41170US	VarseoSmile TriniQ	A1	500 g
41171US	VarseoSmile TriniQ	A2	500 g
41172US	VarseoSmile TriniQ	A3	500 g
41177US	VarseoSmile TriniQ	A3.5	500 g
41173US	VarseoSmile TriniQ	B1	500 g
41178US	VarseoSmile TriniQ	B2	500 g
41174US	VarseoSmile TriniQ	C2	500 g
41179US	VarseoSmile TriniQ	D3	500 g
41175US	VarseoSmile TriniQ	OM1	500 g
41176US	VarseoSmile TriniQ	OM3	500 g
41180US	VarseoSmile TriniQ	A1	250 g
41181US	VarseoSmile TriniQ	A2	250 g
41182US	VarseoSmile TriniQ	A3	250 g
41187US	VarseoSmile TriniQ	A3.5	250 g
41183US	VarseoSmile TriniQ	B1	250 g
41188US	VarseoSmile TriniQ	B2	250 g
41184US	VarseoSmile TriniQ	C2	250 g
41189US	VarseoSmile TriniQ	D3	250 g
41185US	VarseoSmile TriniQ	OM1	250 g
41186US	VarseoSmile TriniQ	OM3	250 g

Shades and shade designations based on VITA Classical A1–D4® shade guide and VITA Toothguide 3D-MASTER®.

Note: The availability of individual product variants may vary from region to region. The latest product information can be found on the BEGO website.

7. Warnings

When processing VarseoSmile TriniQ in liquid form, the following warnings must be observed:

- Read the safety data sheet before use.
- VarseoSmile TriniQ may only be used for the fabrication of definitive or temporary bridges, crowns, inlays, onlays, veneers and denture teeth according to the design requirements described in these Instructions for Use.
- Any deviation from the Instructions for Use may have a negative impact on the chemical and physical properties and biocompatibility of the final product.
- The product contains substances which may cause skin irritation or allergic skin reactions.
- Avoid inhalation or ingestion. May cause irritation to the respiratory tract.
- May damage fertility or the unborn child.
- In case of skin contact with the resin, wash with plenty of water.
- In case of eye contact, rinse cautiously with water for several minutes. Remove contact lenses if present and easily removed. Continue to rinse.
- If swallowed or suspected, seek medical advice and call a poison control center / physician if you feel unwell.

8. Precautionary measures

Wear protective clothing, goggles and nitrile gloves when handling VarseoSmile TriniQ. Further information on the handling of the product can be taken from the safety data sheet available at www.bego.com. It cannot be excluded that in rare cases individual reactions (e.g. incompatibilities or allergic reactions) may occur. In these cases, the user should no longer use VarseoSmile TriniQ.

Store in a well-ventilated place. Keep container tightly closed. Keep under lock and key. Observe the storage and transport conditions. Do not use the resin after the expiration date.

Forward these instructions for use or the relevant information contained therein to other users of the products, as failure to observe them may result in a high risk of unreliable and/or unusable products, which may endanger the safety of the user.

9. MR Safety

VarseoSmile TriniQ is composed of materials that are electrically non-conductive, nonmetallic, and nonmagnetic.

10. Adverse Side-Effects

Adverse effects, reactions or similar occurrences, including those not listed in these instructions for use, should be reported immediately to the manufacturer or local distributor.

11. Transport and storage

VarseoSmile TriniQ resin must be protected from light to prevent spontaneous polymerization. Storage and transport must take place in the original bottle in a temperature range between 4 °C and 35 °C. Resin residues can be stored in the resin tray of the 3D printer with the UV light protection hood of the printer closed or under a light-tight cover of the resin tray. The expiry date up to which the resin can be used is indicated on the label of each original bottle.

Note: Processing the material after the expiration date may lead to unusable results. Fully polymerized printed objects can be stored at room temperature protected from light.

12. Disposal

Completely polymerized material residues (base plate, support structure) can be disposed of in household waste. Unused resin or cleaning solutions with resin residues must be disposed of in accordance with local and national regulations.

13. Processing Instruction

VarseoSmile TriniQ has been developed to produce dental objects using the 3D printing process. To ensure safe and functional results, the following instructions must be observed.

13.1. Equipment Required

Warning: Results that safely meet the indications for use of VarseoSmile TriniQ can only be achieved when using validated 3D printing equipment, cleaning methods and UV polymerization units.

Validated 3D Printers	<ul style="list-style-type: none"> • Varseo XS (405 nm); Firmware 2.6.8.31 and above • Asiga Max UV (385 nm, 405 nm); Firmware 2018-09-03 and above
Validated cleaning method	<ul style="list-style-type: none"> • Ultrasonic bath filled with either: <ul style="list-style-type: none"> • Ethanol 96 % • Isopropanol 99.5 % • InovaPrint Wash*
Validated UV polymerization units	<ul style="list-style-type: none"> • Varseo BEGO Cure (365 nm, 385 nm, 405 nm) • BEGO Ototflash (300 nm – 700 nm)
Further recommended equipment	<ul style="list-style-type: none"> • Stainless steel spatula • Spray bottle with ethanol or isopropanol • Cutting disc or side cutter (for removing the support structure) • Protective equipment • Possibly blasting unit 1.5 bar • Possibly blasting agent 50 µm

An overview of all 3D printers, cleaning agents and UV light-curing units that have been validated and approved for use with VarseoSmile TriniQ can be found in the 3D printing system components compatibility overview at www.begousa.com

13.2. Design requirements

When digitally designing dental objects from VarseoSmile TriniQ, the following specifications must be observed:

Indication	Design	Illustration**	Permanent	Temporary ≤ 12 mos.	Design requirement
Single tooth restorations	Fully anatomical single crown in the anterior region		✓	✓	<ul style="list-style-type: none"> • Minimum wall thickness for posterior and anterior teeth 0.7 mm • Minimum wall thickness for veneers 0.5 mm
	Fully anatomical single crowns in the posterior region		✓	✓	
	Veneer		✓	✓	
	Inlay, onlay		✓	✓	
Bridge constructions	Fully anatomical monolithic bridge constructions for the anterior and/or posterior region, up to 7 units without pontics		✓	✓	<ul style="list-style-type: none"> • Maximum one pontic between two abutment crowns • Bridges in the anterior region: <ul style="list-style-type: none"> • Minimum wall thickness occlusal and circumferential 1.0 mm • Minimum connector cross section 14 mm² • Bridges in the posterior region: <ul style="list-style-type: none"> • Minimum wall thickness occlusal and circumferential 1.0 mm • Minimum connector cross section 16 mm²
	Fully anatomical monolithic bridge constructions for the anterior and/or posterior region with max. one bridge pontic, up to 3 units		✓	✓	
	Fully anatomical monolithic bridge for the anterior and/or posterior region with max. 3 non-adjacent pontics (up to 7-unit bridge constructions)		X	✓	
	Fully anatomical monolithic bridge for the anterior and/or posterior region with a denture base		✓	✓	
Denture teeth	Denture teeth for bonding with a denture base		✓	✓	<ul style="list-style-type: none"> • Minimum wall thickness 1.0 mm

** example

13.3. Nesting

Required software based on the 3D printer used:

- CAM Creator Print version 1.32 or higher
- Asiga Composer version 1.3.3 or higher

The printing parameters required for processing VarseoSmile TriniQ are specified in the specific material file for your 3D printer. To select the correct set of printing parameters, select your 3D printer model and VarseoSmile TriniQ in the respective software. Load the desired STL-File into your software

Further printing parameters:

- Layer thickness 50 µm (preset/default)
- Horizontal alignment of dental objects, occlusal plane to the build platform; Exception for bridges: Orient bridges in an angle of 45° to the build platform
- Use automatic generation of support structures in the respective software to ensure required size, density and spacing of pins

* This symbol is a commercial designation/registered trademark of a company that is not part of the BEGO company group.

13.4. 3D Printing

The ideal processing temperature range of VarseoSmile TriniQ is between 18 °C and 35 °C. Shake the bottle for approx. 2 min before first use to homogenize the resin. When transferring the resin to the resin tray of the validated 3D printer, make sure that the resin is only exposed to daylight for as short a time as possible. Mix the resin in the resin tray when a transparent layer is visible on the surface.

Start the 3D print job according to the 3D printer manufacturer's instructions for use. Remove the printed objects from the build platform according to the 3D printer manufacturer's instructions for use. If necessary, use a stainless-steel spatula for this purpose.

Note: After removal of the printed objects from the 3D printer and removal from the build platform, the subsequent steps of removal of residual resin and post-curing should be performed without prolonged interruption and storage of the printed objects to avoid undesirable material changes.

13.5. Removal of residual resin

Residual resin is removed by cleaning the printed objects in an ultrasonic bath. Ethanol (≥96 %), isopropanol (≥99,5 %) or InovaPrint Wash* can be used as cleaning solution:

1. 3 minutes cleaning in the ultrasonic bath in reusable cleaning solution.
2. 2 minutes cleaning in the ultrasonic bath in fresh cleaning solution.
3. spray the printed objects with ethanol (≥96 %) or isopropanol (≥99,5 %).
4. Carefully dry the cleaned objects with compressed air.
5. Remaining shiny spots on the printed objects indicate remaining residual resin. Remove this by spraying the printed objects with ethanol (≥96 %) or isopropanol (≥99,5 %). Residual resin can also be removed with a brush soaked in ethanol (≥96 %) or isopropanol (≥99,5 %).
6. If necessary, dry the cleaned objects again with compressed air

Warning: Never fill ethanol or isopropanol directly into the ultrasonic bath, but always place it in the recommended container (REF 19621) in the ultrasonic bath filled with water. An explosion-proof ultrasonic bath must be used.

Note: The specified total cleaning time should not be exceeded, otherwise the printed objects may be damaged.

Note: Especially when using ethanol (≥96 %) or isopropanol (≥99,5 %) as a cleaning solution, a whitish layer may form on the surface of the dental objects. This is a deposit of the ceramic filling of VarseoSmile TriniQ. The deposit is harmless and can be removed by polishing after post-curing of the printed objects.

13.6. Post-curing

Remove the support structures with a cutting disk or side cutter before post-curing. Take care not to deform the printed objects.

Post-cure the printed objects as shown in following table:

Curing unit	Curing parameters	Curing process
BEGO Varseo Cure	Choose "VarseoSmile TriniQ" (2 x 2:30 min)	Turn object between cycles
BEGO Otoflash	2 x 2000 Flashes	Turn object between cycles

An overview of other UV light-curing units that have been validated and approved for post-curing dental restorations made from VarseoSmile TriniQ can be found in the 3D printing system components compatibility overview at www.begousa.com

Note: The duration of post-curing specified in each case applies only to regularly maintained UV curing units.

Note: After post-curing, a slight yellowish discoloration of the printed objects may appear. This discoloration is normal and will disappear after a short time.

Warning: If there is an unintended interruption of the post-curing process, do not use the printed object until it has been polymerized under a full post-curing cycle.

13.7. Finishing and further processing of printed objects

Post-cured objects made of VarseoSmile TriniQ can be polished with pumice stone and polishing paste. Take care not to deform the printed objects when polishing.

Post-cured objects made of VarseoSmile TriniQ can be characterized or finished and supplemented with composite stains as well as direct and indirect composites. The valid instructions for use of the manufacturer of the composite stains or composite veneering material must be observed.

Cementation of bridges, crowns, inlays, onlays and veneers

- Permanent cementation of objects made of VarseoSmile TriniQ can be performed with fully adhesive or self-adhesive luting composites.
- Temporary cementation of objects made of VarseoSmile TriniQ can be performed with temporary luting cements.
- The cementation of objects made of VarseoSmile TriniQ on metal frameworks can be performed with adhesive luting or veneering composites.

Note: In any case, the instructions for use of the respective composite or cement manufacturer must be observed.

Note: Glass ionomer cements are not suitable for cementing objects made of VarseoSmile TriniQ.

Fixation of denture teeth

Denture teeth made of VarseoSmile TriniQ can be bonded to cold-polymerized or 3D-printed denture base resins. For this purpose, the respective valid instructions of the denture base material manufacturer for bonding denture teeth to the denture base must be observed. If the denture base manufacturer does not provide instructions for bonding the denture teeth to the denture base, the denture base material can be used to bond the denture teeth to the denture base in the case of cold-cured denture bases and the liquid VarseoSmile TriniQ resin can be used to bond the denture teeth to the denture base in the case of 3D-printed denture bases.

Cleaning of post-cured objects

Fully post-cured objects made of VarseoSmile TriniQ can be cleaned by evaporation or in an immersion bath with cleaning solution.

14. Symbols used on labels

Reference number and graphic	Title	Description	Reference
5.1.1 	Manufacturer	Indicates the manufacturer of the medical device	ISO 15223-1:2021 Medical devices • Symbols to be used with information to be supplied by the manufacturer
5.1.3 	Date of manufacture	Indicates the date when the medical device was manufactured	ISO 15223-1:2021 Medical devices • Symbols to be used with information to be supplied by the manufacturer
5.1.4 	Use-by date	Indicates the date after which the medical device is not to be used	ISO 15223-1:2021 Medical devices • Symbols to be used with information to be supplied by the manufacturer
5.1.5 	Batch code	Indicates the manufacturer's batch code so that the batch or lot can be identified	ISO 15223-1:2021 Medical devices • Symbols to be used with information to be supplied by the manufacturer
5.1.6 	Catalogue number	Indicates the manufacturer's catalogue number so that the medical device can be identified	ISO 15223-1:2021 Medical devices • Symbols to be used with information to be supplied by the manufacturer
5.3.2 	Keep away from sunlight	Indicates a medical device that needs protection from light sources	ISO 15223-1:2021 Medical devices • Symbols to be used with information to be supplied by the manufacturer
5.3.7 	Temperature limit	Indicates the temperature limits to which the medical device can be safely exposed	ISO 15223-1:2021 Medical devices • Symbols to be used with information to be supplied by the manufacturer
5.4.3 	Consult instructions for use or consult electronic instructions for use	Indicates the need for the user to consult the instructions for use	ISO 15223-1:2021 Medical devices • Symbols to be used with information to be supplied by the manufacturer
5.4.4 	Caution	Indicates that caution is necessary when operating the device or control close to where the symbol is placed, or that the current situation needs operator awareness or operator action in order to avoid undesirable consequences	ISO 15223-1:2021 Medical devices • Symbols to be used in the information provided by the manufacturer
	For professional only	Caution: Federal law restricts this device to sale by or on the order of a dentist	FDA Use of Symbols in Labeling final rule
	MR safe	Indicates that the device is safe to use in magnetic resonance environment	ASTM F2503-20

* This symbol is a commercial designation/registered trademark of a company that is not part of the BEGO company group.

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